




**Molecular Epidemiology, Inc.**  
**Sample Analysis Report**

dba   
 Laboratories & Consulting Group

F09/MB450/09

<b>Customer ID Number</b>	XXXX-XXXX	<b>Sample Received</b>	00/00/0000	<b>MEI Report #</b>	XXXXXXXXXX-XX	<b>Date Reported</b>	00/00/0000
---------------------------	-----------	------------------------	------------	---------------------	---------------	----------------------	------------

<b>Customer Sample #</b>	XXXXX	<b>Job ID</b>	XXXXX-XXXX-XXXX
<b>Sample Type</b>	Agar plate	<b>Requested Service(s)</b>	Phenotypic and Genetic Identification
<b>Sample Source</b>	Contaminant	<b>Test Method(s)</b>	MB450, MB457, MB473
<b>Customer Comments</b>	None		
<b>MEI Sample #</b>	XXXXX-X		

Test	Test Results	
<b>Gram Stain Reaction</b>	Negative	
<b>Microscopic</b>	Rods	
<b>Oxidase</b>	Negative	
<b>Catalase</b>	Positive	

<b>Genetic ID</b> Comparisons to genetically similar microorganisms			<b>Notes:</b> Supplemental phenotypic data support <i>Burkholderia cepacia</i> complex. <i>B.cepacia</i> , <i>B. anthina</i> , <i>B. cenocepacia</i> , and <i>B. multivorans</i> belong to this complex.
<b>Genetic Distance</b>	<b>Genus</b>	<b>Species</b>	
0.0001	<i>Burkholderia</i>	<i>cepacia</i>	
0.0002	<i>Burkholderia</i>	<i>cenocepacia</i>	
0.0005	<i>Burkholderia</i>	<i>anthina</i>	
0.0022	<i>Burkholderia</i>	<i>multivorans</i>	

<b>Deviations</b>	None
-------------------	------

<b>Microbial ID Conclusion</b>	<i>Burkholderia cepacia</i> complex
--------------------------------	-------------------------------------

<b>Microbiology Review</b>	_____	<b>Date:</b> 00/00/0000
<b>Genetics Review</b>	_____	<b>Date:</b> 00/00/0000
<b>QC Approval</b>	_____	<b>Date:</b> 00/00/0000
<b>QA Approval</b>	_____	<b>Date:</b> 00/00/0000

<b>Reported By:</b>	<b>Molecular Epidemiology, Inc.</b>
Phone: 206-306-8882	15300 Bothell Way NE,
Fax: 206-306-8883	Lake Forest Park, WA 98155
e-Mail: <a href="mailto:sidhuj@molecularepi.com">sidhuj@molecularepi.com</a>	<a href="http://www.molecularepi.com">www.molecularepi.com</a>

This laboratory report relates only to the portion of the sample which was tested in this report. Interpretation of these results is the sole responsibility of the Customer.

TRADE SECRET - CONTAINS CONFIDENTIAL COMMERCIAL INFORMATION

This document contains confidential commercial information pursuant to 5 U.S.C sec. 552(b)(4). Copyright 2007 © by MEI